

Dry Suction/Dry Seal Chest Drain

Product Code: ST382



Advantages of Vide's Dry Suction/Dry Seal Versus Vide's Water Seal Chest Drain

- ✓ **Intuitive Sliding Suction Control panel** - no more need for finding and filling sterile water!
- ✓ **Knock-over nozzle** to ensure Function is retained even after unit is accidentally knocked over
- ✓ **7-Scale Air Leak Monitor** for more precise air leak quantification
- ✓ **8-Discrete Suction Levels** with Suction Control Indicator
 - 5cmH₂O, -10cmH₂O, -15cmH₂O, -20cmH₂O
 - 25cmH₂O, -30cmH₂O, -35cmH₂O & -40cmH₂O
- ✓ **Y-Shape Sampling Port** with self-sampling needless feature, allowing access to fresh blood samples
- ✓ **Low Positive Pressure Relief Valve** release automatically at 2cmH₂O
- ✓ **Negative Pressure Relief Valve** release automatically during prolonged episode of negative intrapleural pressure
- ✓ **Suction and Negative Pressure Indicators** allow easy visualization of the status of the unit
- ✓ Vide's **Calibrated** Dry Suction Chest Drain can maintain consistent suction control at varying negative pressure or flow rate from the source



ALLEVA MEDICAL CO.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong
Tel: +852 3166 7200 Fax: +852 2355 7663 www.allevamedical.com

